

# PolyBLOCK Plate


## For Any Lab Based Screening Offering Even, Controlled And Reproducible Mixing

The PolyBLOCK Plate is designed for use with all standard SBS footprint plates and features two independently controlled heated (up to 190 °C) positions with magnetic stirring.


HEL has extensively checked the magnetic stirring on this system and it provides effective magnetic agitation in all positions on a standard plate.

Magnetic agitation is superior to mixing achieved by vortexing and avoids the problems associated with the strong magnetic fields used for tumble stirring.

The PolyBLOCK Plate makes an excellent addition to any lab doing plate based screening. It gives even, controlled and reproducible mixing and heating on a fully automated platform.



**Easy To Use**



**Magnetic Agitation**



# About HEL

HEL is an international company that specialises in chemical reactors, bioreactors and related data/logging tools for process R&D in the pharmaceutical, fine chemical, biotechnology and petrochemical industries.

Established in 1987 and with the clients worldwide our key strengths are:

Knowledgeable staff – highly qualified and experienced chemical engineers and chemists

Quality – underpinned by ISO9001 certification for over 16 years

Service – choice of service contracts backed by an established culture of unmatched client support

Range of products – both of-the-shelf and custom designs, manual and fully automated controls, low and high pressure/temperature applications, single and parallel/multi-vessel products



HEL Ltd // UK

HEL Inc // US

HEL China // Beijing

HEL Germany // Hamburg

HEL India // Thane

HEL Italia // Milan

HEL Russia // Moscow

e: [marketing@helgroup.com](mailto:marketing@helgroup.com)

e: [marketing@hel-inc.com](mailto:marketing@hel-inc.com)

e: [info@helchina.com](mailto:info@helchina.com)

e: [dach@helgroup.com](mailto:dach@helgroup.com)

e: [info@helindia.com](mailto:info@helindia.com)

e: [helitalia@helgroup.com](mailto:helitalia@helgroup.com)

e: [marketing@helgroup.com](mailto:marketing@helgroup.com)

t: +44 20 8736 0640

t: +1 609 912 1551

t: +86 10 8210 1033

t: +44 20 8736 0640

t: +91 22 4153 3317

t: +44 20 8736 0640

t: +44 20 8736 0640

**HEL Limited, 9 - 10 Capital Business Park, Manor Way, Borehamwood, Herts  
WD6 1GW, UK**



better chemistry – faster

[www.helgroup.com](http://www.helgroup.com)